









# Table of Contents

|                    |    |
|--------------------|----|
| Indications        | 2  |
| Design Features    | 3  |
| Templating         | 4  |
| Patient Prep       | 5  |
| Surgical Technique | 6  |
| Catalog            | 18 |

## Setting the Standard Since 1956

For more than 45 years, the Compression Hip Screw from Smith & Nephew has set the standard for treating femoral fractures.

CHS applications include intertrochanteric, subtrochanteric, subcapital, and intracapsular neck fractures. A wide variety of implant sizes extends the CHS versatility, and the unsurpassed instrumentation is the result of years of research and clinical testing.

The CHS is a complete system, offering surgeons many unique design features:

- Keyed and keyless plates
- Cold-forged plates for improved fatigue strength
- Deeper threads in the lag screw for increased bone purchase
- I-beam construct, improving fatigue strength in the lag screw
- Adult, intermediate, and pediatric versions

Maximize your options in surgery and rely on the proven performance and consistent results of the Compression Hip Screw from Smith & Nephew.



# Indications

The Smith & Nephew Compression Hip Screw is intended to treat fractures of the proximal femur.

This includes:

- Intertrochanteric Fractures
- Subtrochanteric Fractures
- Intracapsular Neck Fractures
- Subcapital Fractures





# Design Features

## CLASSIC Keyed and AMBI® Keyed/Keyless Plates

**AMBI Plates:** Barrel design is keyless but can be converted to keyed with the insertion of a small keying clip.


**CLASSIC Plates:** Barrel design is keyed only.

Angles: 130° to 150° in 5° increments

Lengths: 60mm to 300mm, 2 to 14 slots

Barrel Lengths: Standard plate barrels are 38.1mm long

Selected sizes are available with a shorter 25.9mm barrel



### AMBI CLASSIC Lag Screws

18 Lengths: 55mm–140mm

Nonself-tapping for cancellous bone

Scratch Resistant Surface (SRS)

### Standard Lag Screws

Thread Diameter: 12.7mm

Root Diameter: 9.0mm

### Super Lag Screws

Thread Diameter: 14.3mm

Root Diameter: 9.0mm

Medial fragments can be captured by inserting a 6.5mm cannulated or 6.5mm cancellous screw through the proximal slot. 6.5mm cannulated screws are available from 25mm to 120mm in 5mm increments.

Compression slots can be used to gain up to 2mm of compression along the shaft of the femur.

The distal slot in plates with 8 or more holes can be used to gain up to 5.5mm of shaft compression.

### AMBI CLASSIC Compression Screws

Lengths – 19.0mm and 28.5mm

Hex Diameter – 3.5mm

### 4.5mm Self-tapping Cortical Bone

Screws: 25 Lengths – 16mm–64mm

Major Thread Diameter – 4.5mm

Minor Thread Diameter – 3.2mm


Hex Diameter – 3.5mm



# Templating

Smith & Nephew offers preoperative X-ray templates (7118-0454) to ensure selection of the appropriate implant for each fracture. Templating helps determine the proper lag screw length as well as the correct neck angle and size of the plate.

Templates are available at a magnification of 117% of actual size. An anteroposterior pelvic X-ray and a cross table lateral view of the hip are necessary when templating. For best results, templates should be used on X-rays that feature the extremity placed in slight internal rotation with gentle longitudinal traction applied.






# Patient Prep

## Preoperative Planning

Adequate preoperative assessment of the fracture requires a thorough history, careful physical exam, and adequate radiographic studies. Ambulatory status as well as previous lower extremity fractures and surgeries must be known. Obtain appropriate radiographic images. Note any skin compromise.



## Patient Positioning

Patient is placed supine on a fracture table with the unaffected leg flexed at the hip and knee, and then abducted and slightly internally rotated.



# Surgical Technique

## Guide Pin Insertion

Before inserting the guide pin, satisfactory fracture reduction using an open or closed technique should be accomplished. Using the multiple angle guide or the adjustable angle guide, aim the 3.2mm guide pin toward the apex of the femoral head, ensuring that it is parallel to and in the center of the femoral neck. Advance the guide pin to within 10mm of the joint line to prevent joint penetration.




Multiple Angle Guide



Adjustable Angle Guide



Quick Connect Adaptor



Tip Threaded Guide Pin

Percutaneous Direct Measuring Gauge





Guide Pin Placement Instrument

When the guide pin is adequately positioned in the femoral head, use the percutaneous direct measuring gauge to determine the appropriate lag screw length and reaming distance.

**Note:** This instrument directly measures the guide pin.

The guide pin placement instrument will allow the placement of a parallel guide pin 13mm proximal to the primary guide pin. This will maintain reduction of unstable fractures during reaming.








## Reaming


Set the power combination reamer to the lag screw length indicated by the measuring gauge. Ream for lag screw using the combination reamer. Insert the reamer over the guide pin. To minimize the occurrence of guide pin pullout, it is important to avoid reaming over the threaded portion of the guide pin. This can be achieved in two ways.



Short barrel notch indicator to be used when reaming short barrel plates.



Power Combination Reamer



Option 1: Once the guide pin is inserted and measured, advance it an additional 5mm into the subchondral bone and ream according to the exact lag screw length measurement. Choose a lag screw that matches the length measurement.

Option 2: Insert the guide pin into the subchondral bone, measure, and set the reamer 5mm shorter than the length measured. Choose a lag screw that matches the length that was reamed.



Quick-Connect T-Handle

## Tapping

Tapping is indicated in younger patients or abnormally dense bone. It is also indicated to avoid excessive torque during insertion of the lag screw. To tap, attach the quick-connect T-handle to the lag screw tap and set it to the appropriate lag screw length. Tap until the advancing portion of the positive stop rests against the cortex guide.



Lag Screw Tap



## Lag Screw Length Selection

Lag screw length will determine the amount of compression you can achieve with the compression screw intraoperatively. If the lag screw selected is exactly the same length of the measurement taken from the direct measuring gauge, this will provide 5mm of compression. If more compression is needed, a 5mm shorter lag screw will permit 5mm of additional compression.



# CLASSIC Wrench and Plate Assembly Technique

To Be Used with CLASSIC Plates Only.



Slide the CLASSIC insertion wrench through the barrel of the plate and attach the lag screw.

**Do Not Use Wrench as a Lever.**

## Assembled



Slip the AMBI CLASSIC centering sleeve onto the insertion wrench.




AMBI CLASSIC  
Centering Sleeve

## Lag Screw Insertion

Advance the lag screw to the predetermined depth and verify using image intensification. As a guide, the rings are aligned on the AMBI CLASSIC centering sleeve and CLASSIC insertion wrench and these can be used to aid in determining lag screw seating. These marks are calibrated for 135° plates and 150° plates.

When the lag screw is inserted to the desired depth, the wrench handle must be perpendicular to the axis of the femoral shaft for proper keying to the plate barrel. Remove the AMBI CLASSIC centering sleeve and advance the plate onto the lag screw shaft.



CLASSIC  
Insertion Wrench



## Final Seating of the Plate

Use a plate tamper to fully seat the plate.

Unscrew the retaining rod from the lag screw and remove the insertion wrench from the back of the lag screw. Remove the guide pin.



Plate Tamper

## Plate Tamper

Unscrew the retaining rod from the lag screw and remove the insertion wrench. Slide the cannulated plate tamper over the guide pin into the barrel of the plate to fully seat the plate. Remove the guide pin once the plate is fully seated.



Cannulated Plate Tamper



# AMBI Wrench, Plate and Clip Assembly Technique

To Be Used with AMBI Plates Only.



Assemble the AMBI clip, plate and lag screw onto the AMBI insertion wrench. The AMBI clip may be omitted for a keyless system.

**Do Not Use Wrench as a Lever.**

Slip the AMBI CLASSIC centering sleeve onto the insertion wrench.

## Assembled




AMBI CLASSIC  
Centering Sleeve

## Lag Screw Insertion


Advance the lag screw to the predetermined depth and verify using image intensification. As a guide, when using the 135° plate, the lag screw should be advanced until the rings on the AMBI insertion wrench are aligned with the 135° marks on the AMBI CLASSIC centering sleeve. This process remains the same for the 150° plate.

When the lag screw is inserted to the desired depth, the wrench handle must be perpendicular to the axis of the femoral shaft for proper keying to the plate barrel. Remove the AMBI CLASSIC centering sleeve and advance the plate onto the lag screw shaft.



AMBI  
Insertion Wrench





## Clip Insertion

For Use with AMBI Plates Only.

Make sure the wrench is perpendicular to the femur and be sure the lines on the plate barrel and AMBI insertion wrench are aligned. Press the AMBI clip inserter down the shaft of the AMBI insertion wrench. Use firm finger pressure until the clip snaps in place.



AMBI Clip



## Final Seating of the Plate

Use a plate tamper to fully seat the plate.

Unscrew the retaining rod from the lag screw and remove the insertion wrench from the back of the lag screw. Remove the guide pin.



## Plate Tamper


Unscrew the retaining rod from the lag screw and remove the insertion wrench. Slide the cannulated plate tamper over the guide pin into the barrel of the plate to fully seat the plate. Remove the guide pin once the plate is fully seated.





# Cannulated Barrel Guide and Insertion/Removal Wrench

To be Used with both CLASSIC and AMBI Plates.

## Wrench Assembly



Attach the cannulated barrel guide to the lag screw. Slide the centering sleeve onto the insertion/removal wrench. Place the entire assembly over the guide pin and advance. **Do Not Use Wrench as a Lever.**



Centering Sleeve



Insertion/Removal Wrench

## Lag Screw Insertion


Advance the lag screw to the desired length and verify using image intensification. As a guide, the rings are aligned on the centering sleeve and insertion/removal wrench and these can be used to aid in determining lag screw seating. These marks are calibrated for 135° plates and 150° plates. Once final depth is obtained, the handle of the insertion/removal wrench must be perpendicular to the axis of the femoral shaft for proper alignment.



## Final Seating of the Plate

Use a plate tamper to fully seat the plate.

Unscrew the retaining rod from the lag screw and remove the insertion wrench from the back of the lag screw. Remove the guide pin.



## Plate Tamper

Unscrew the retaining rod from the lag screw and remove the insertion wrench. Slide the cannulated plate tamper over the guide pin into the barrel of the plate to fully seat the plate. Remove the guide pin once the plate is fully seated.








Plate Clamp

## Attaching the Plate

Use the plate clamp to secure the plate to the shaft.

Use the 3.5mm twist drill through the combination drill guide to drill the bone plate screw holes. Determine appropriate screw length using the bone screw length gauge.




3.5 Twist Drill



Bone Screw Length Gauge

Secure the plate to the bone using 4.5mm self-tapping bone screws.

Insert the remaining screws into the shaft of the femur. Remove the plate clamp.







## Inserting the Compression Screw

Two compression screws are available:


19mm Length – allows 5mm of compression intraoperatively.

28.5mm Length – allows up to 10mm of compression when used in combination with a lag screw that is 5mm shorter than the lag screw measurement.

**Note:** Be sure to always remove the 28.5mm length compression screw and replace with the 19mm length to ensure proper lag screw and plate interface.



Self-Holding  
Hex  
Screwdriver



Final tightening of the compression screw can be achieved with a hex screwdriver once all the 4.5mm screws have been inserted into the plate and traction has been released. Fracture compression is accomplished by fully seating the lag screw with the hex screwdriver.

**Note:** Be sure to take into consideration bone quality and fracture stability while inserting the compression screw to prevent over tightening. Overtightening could cause the lag screw to pull out of the femoral head.



Hex Screwdriver




## Special Features of the Plate

The oval autocompression holes of the plate will allow for up to 2mm of fracture line compression for subtrochanteric fractures or osteotomies.

Place the eccentric gold end of the combination drill guide in the first compression slot distal to the fracture with the arrow pointing towards the fracture. Drill with the 3.5mm Twist drill


Place a 4.5mm self-tapping cortical screw in the slot.




As the screw is seated, it abuts the inclined distal aspect of the slot, forcing the plate and the attached proximal fragment distally until resisted by fracture compression.

For an additional 1mm of compression, repeat this step in the compression slot distal to the first one. Slightly loosen the first eccentrically placed screw after the second screw abuts the slot, but before it is fully seated to allow the additional compression.

Following seating of the second screw, retighten the first screw.







## Capturing the Lesser Trochanter and Posterior-Medial Fragments

The most proximal slot in the plate allows for the insertion of a 6.5mm cancellous or cannulated screw.

These screws can be used to capture the lesser trochanter or a large posterior-medial fragment.

The slot allows up to 45° of proximal and 26° of distal angulation in the coronal plane, and 14° of anterior or posterior angulation in the sagittal plane.



2.4mm Pin Guide

## Inserting the Cannulated Screw

Snap the 2.4mm pin guide into the combination drill guide. Place the pin guide into the proximal slot and insert a 2.4mm guide pin.

2.4mm Guide Pin

## Inserting the 6.5mm Cannulated Screw

**Note:** You will need the instruments from the cannulated screw set to insert the 6.5mm cannulated screw.




Lag Screw Trephine

## Removing the Hip Screw

Remove the compression screw. Next remove the cortical screws from the plate. Lift the plate from the femoral shaft.

Use the lag screw trephine to remove any tissue or bone that may interfere with the lag screw removal.

Connect the insertion/removal wrench to the base of the lag rod via the retaining rod. Using a counterclockwise motion, remove the lag screw.



Retaining Rod for Insertion/Removal Wrench



Insertion/Removal Wrench





CLASSIC  
Compression Hip Screw  
Standard Barrel Plates

Barrel Length: 1.5" (38.1mm)

| Cat. No. | Description | Length | Angle |
|----------|-------------|--------|-------|
| 12-4120  | 2 Slot      | 60mm   | 130°  |
| 12-4121  | 2 Slot      | 60mm   | 135°  |
| 12-4122  | 2 Slot      | 60mm   | 140°  |
| 12-4123  | 2 Slot      | 60mm   | 145°  |
| 12-4124  | 2 Slot      | 60mm   | 150°  |
| 12-4125  | 3 Slot      | 80mm   | 130°  |
| 12-4126  | 3 Slot      | 80mm   | 135°  |
| 12-4127  | 3 Slot      | 80mm   | 140°  |
| 12-4128  | 3 Slot      | 80mm   | 145°  |
| 12-4129  | 3 Slot      | 80mm   | 150°  |
| 12-4130  | 4 Slot      | 100mm  | 130°  |
| 12-4131  | 4 Slot      | 100mm  | 135°  |
| 12-4132  | 4 Slot      | 100mm  | 140°  |
| 12-4133  | 4 Slot      | 100mm  | 145°  |
| 12-4134  | 4 Slot      | 100mm  | 150°  |
| 12-4135  | 5 Slot      | 120mm  | 130°  |
| 12-4136  | 5 Slot      | 120mm  | 135°  |
| 12-4137  | 5 Slot      | 120mm  | 140°  |
| 12-4138  | 5 Slot      | 120mm  | 145°  |
| 12-4139  | 5 Slot      | 120mm  | 150°  |
| 12-4140  | 6 Slot      | 140mm  | 130°  |
| 12-4141  | 6 Slot      | 140mm  | 135°  |
| 12-4142  | 6 Slot      | 140mm  | 140°  |
| 12-4143  | 6 Slot      | 140mm  | 145°  |
| 12-4144  | 6 Slot      | 140mm  | 150°  |
| 12-4145  | 8 Slot      | 180mm  | 130°  |
| 12-4146  | 8 Slot      | 180mm  | 135°  |
| 12-4147  | 8 Slot      | 180mm  | 140°  |
| 12-4148  | 8 Slot      | 180mm  | 145°  |
| 12-4149  | 8 Slot      | 180mm  | 150°  |
| 12-4150  | 10 Slot     | 220mm  | 130°  |
| 12-4151  | 10 Slot     | 220mm  | 135°  |
| 12-4152  | 10 Slot     | 220mm  | 140°  |
| 12-4153  | 10 Slot     | 220mm  | 145°  |
| 12-4154  | 10 Slot     | 220mm  | 150°  |
| 12-4156  | 12 Slot     | 260mm  | 135°  |
| 12-4158  | 12 Slot     | 260mm  | 145°  |
| 12-4161  | 14 Slot     | 300mm  | 135°  |
| 12-4163  | 14 Slot     | 300mm  | 145°  |



CLASSIC  
Compression Hip Screw  
Short Barrel Plates

Barrel Length: 1.0" (25.9mm)

| Cat. No. | Description | Length | Angle |
|----------|-------------|--------|-------|
| 12-4176  | 4 Slot      | 100mm  | 130°  |
| 12-4177  | 4 Slot      | 100mm  | 135°  |
| 12-4178  | 4 Slot      | 100mm  | 140°  |
| 12-4179  | 4 Slot      | 100mm  | 145°  |
| 12-4180  | 4 Slot      | 100mm  | 150°  |
| 12-4165  | 5 Slot      | 120mm  | 130°  |
| 12-4166  | 5 Slot      | 120mm  | 135°  |
| 12-4167  | 5 Slot      | 120mm  | 140°  |
| 12-4168  | 5 Slot      | 120mm  | 145°  |
| 12-4169  | 5 Slot      | 120mm  | 150°  |



## AMBI

### Compression Hip Screw Standard Barrel Plates

Barrel Length: 1.5" (38.1mm)

| Cat. No. | Description | Length | Angle |
|----------|-------------|--------|-------|
| 12-1120  | 2 Slot      | 60mm   | 130°  |
| 12-1121  | 2 Slot      | 60mm   | 135°  |
| 12-1122  | 2 Slot      | 60mm   | 140°  |
| 12-1123  | 2 Slot      | 60mm   | 145°  |
| 12-1124  | 2 Slot      | 60mm   | 150°  |
| 12-1125  | 3 Slot      | 80mm   | 130°  |
| 12-1126  | 3 Slot      | 80mm   | 135°  |
| 12-1127  | 3 Slot      | 80mm   | 140°  |
| 12-1128  | 3 Slot      | 80mm   | 145°  |
| 12-1129  | 3 Slot      | 80mm   | 150°  |
| 12-1130  | 4 Slot      | 100mm  | 130°  |
| 12-1131  | 4 Slot      | 100mm  | 135°  |
| 12-1132  | 4 Slot      | 100mm  | 140°  |
| 12-1133  | 4 Slot      | 100mm  | 145°  |
| 12-1134  | 4 Slot      | 100mm  | 150°  |
| 12-1135  | 5 Slot      | 120mm  | 130°  |
| 12-1136  | 5 Slot      | 120mm  | 135°  |
| 12-1137  | 5 Slot      | 120mm  | 140°  |
| 12-1138  | 5 Slot      | 120mm  | 145°  |
| 12-1139  | 5 Slot      | 120mm  | 150°  |
| 12-1140  | 6 Slot      | 140mm  | 130°  |
| 12-1141  | 6 Slot      | 140mm  | 135°  |
| 12-1142  | 6 Slot      | 140mm  | 140°  |
| 12-1143  | 6 Slot      | 140mm  | 145°  |
| 12-1144  | 6 Slot      | 140mm  | 150°  |
| 12-1145  | 8 Slot      | 180mm  | 130°  |
| 12-1146  | 8 Slot      | 180mm  | 135°  |
| 12-1147  | 8 Slot      | 180mm  | 140°  |
| 12-1148  | 8 Slot      | 180mm  | 145°  |
| 12-1149  | 8 Slot      | 180mm  | 150°  |
| 12-1150  | 10 Slot     | 220mm  | 130°  |
| 12-1151  | 10 Slot     | 220mm  | 135°  |
| 12-1152  | 10 Slot     | 220mm  | 140°  |
| 12-1153  | 10 Slot     | 220mm  | 145°  |
| 12-1154  | 10 Slot     | 220mm  | 150°  |
| 12-1156  | 12 Slot     | 260mm  | 135°  |
| 12-1158  | 12 Slot     | 260mm  | 145°  |
| 12-1161  | 14 Slot     | 300mm  | 135°  |
| 12-1163  | 14 Slot     | 300mm  | 145°  |







## AMBI Compression Hip Screw Short Barrel Plates

Barrel Length: 1.0" (25.4mm)

| Cat. No. | Description | Length | Angle |
|----------|-------------|--------|-------|
| 12-1198  | 4 Slot      | 100mm  | 130°  |
| 12-1199  | 4 Slot      | 100mm  | 135°  |
| 12-1200  | 4 Slot      | 100mm  | 140°  |
| 12-1201  | 4 Slot      | 100mm  | 145°  |
| 12-1202  | 4 Slot      | 100mm  | 150°  |
| 12-1165  | 5 Slot      | 120mm  | 130°  |
| 12-1166  | 5 Slot      | 120mm  | 135°  |
| 12-1167  | 5 Slot      | 120mm  | 140°  |
| 12-1168  | 5 Slot      | 120mm  | 145°  |
| 12-1169  | 5 Slot      | 120mm  | 150°  |



## AMBI CLASSIC/IMHS Standard Lag Screws

Thread Diameter: 1/2" (12.7mm)

Thread Length: 21.0mm

Root Diameter: 9.0mm

| Cat. No. | Length | Cat. No. | Length |
|----------|--------|----------|--------|
| 12-1100  | 55mm   | 12-1109  | 100mm  |
| 12-1101  | 60mm   | 12-1110  | 105mm  |
| 12-1102  | 65mm   | 12-1111  | 110mm  |
| 12-1103  | 70mm   | 12-1112  | 115mm  |
| 12-1104  | 75mm   | 12-1113  | 120mm  |
| 12-1105  | 80mm   | 12-1114  | 125mm  |
| 12-1106  | 85mm   | 12-1176  | 130mm  |
| 12-1107  | 90mm   | 12-1177  | 135mm  |
| 12-1108  | 95mm   | 12-1178  | 140mm  |

NOTE: Do not use AMBI/CLASSIC 55, 60, or 65mm lag screws with IMHS. These sizes are too short to work effectively with this device.



## AMBI Clip Cat. No. 12-1115



## AMBI CLASSIC Super Lag Screws

Thread Diameter: 9/16" (14.3mm)

Thread Length: 21.0mm

Root Diameter: 9.0mm

| Cat. No. | Length | Cat. No. | Length |
|----------|--------|----------|--------|
| 12-1180  | 55mm   | 12-1189  | 100mm  |
| 12-1181  | 60mm   | 12-1190  | 105mm  |
| 12-1182  | 65mm   | 12-1191  | 110mm  |
| 12-1183  | 70mm   | 12-1192  | 115mm  |
| 12-1184  | 75mm   | 12-1193  | 120mm  |
| 12-1185  | 80mm   | 12-1194  | 125mm  |
| 12-1186  | 85mm   | 12-1195  | 130mm  |
| 12-1187  | 90mm   | 12-1196  | 135mm  |
| 12-1188  | 95mm   | 12-1197  | 140mm  |

NOTE: Not compatible with IMHS nail.



## 6.5mm Cannulated Screws

| Cat. No. | Length | Cat. No. | Length |
|----------|--------|----------|--------|
| 12-1625  | 25mm   | 12-1635  | 75mm   |
| 12-1626  | 30mm   | 12-1636  | 80mm   |
| 12-1627  | 35mm   | 12-1637  | 85mm   |
| 12-1628  | 40mm   | 12-1638  | 90mm   |
| 12-1629  | 45mm   | 12-1639  | 95mm   |
| 12-1630  | 50mm   | 12-1640  | 100mm  |
| 12-1631  | 55mm   | 12-1641  | 105mm  |
| 12-1632  | 60mm   | 12-1642  | 110mm  |
| 12-1633  | 65mm   | 12-1643  | 115mm  |
| 12-1634  | 70mm   | 12-1644  | 120mm  |



## Compression Screws

| Cat. No. | Length | Cat. No. | Length |
|----------|--------|----------|--------|
| 12-1116  | 19mm   | 12-1117  | 28.5mm |



## 4.5mm Self-Tapping Cortical Bone Screws

Head Diameter: 8.0mm

Major Thread Diameter: 4.5mm

Root Diameter: 3.2mm

| Cat. No.  | Length | Cat. No.  | Length |
|-----------|--------|-----------|--------|
| 7112-9216 | 16mm   | 7112-9242 | 42mm   |
| 7112-9218 | 18mm   | 7112-9244 | 44mm   |
| 7112-9220 | 20mm   | 7112-9246 | 46mm   |
| 7112-9222 | 22mm   | 7112-9248 | 48mm   |
| 7112-9224 | 24mm   | 7112-9250 | 50mm   |
| 7112-9226 | 26mm   | 7112-9252 | 52mm   |
| 7112-9228 | 28mm   | 7112-9254 | 54mm   |
| 7112-9230 | 30mm   | 7112-9256 | 56mm   |
| 7112-9232 | 32mm   | 7112-9258 | 58mm   |
| 7112-9234 | 34mm   | 7112-9260 | 60mm   |
| 7112-9236 | 36mm   | 7112-9262 | 62mm   |
| 7112-9238 | 38mm   | 7112-9264 | 64mm   |
| 7112-9240 | 40mm   |           |        |







### Tip Threaded Guide Pin, 3.2mm

| Cat. No.  | Description                |
|-----------|----------------------------|
| 7111-0056 | Sterile Package, Single    |
| 11-0016   | Nonsterile Package, Single |



Quick-Connect Adaptor  
(for use with Hall, Jacobs, or Stryker power)  
Cat. No. 7111-0020



Quick-Connect Adaptor  
(for use with Synthes power)  
Cat. No. 7111-0022



Perforation Drill  
Cat. No. 11-0021  
NOTE: Used to make opening in lateral cortex guide pin.



Multiple Angle Guide  
Cat. No. 7111-0028



Quick-Connect T-Handle  
Cat. No. 7111-5045



Adjustable Angle Guide  
Cat. No. 11-0925



Percutaneous Direct Measuring Gauge  
Cat. No. 11-0026



## Guide Pin Placement Instrument

Cat. No. 11-5036



## Power Combination Reamer

Cat. No. 11-0023



## Trial Handle

Cat. No. 11-0047



## Trial Plates

| Cat. No. | Angle |
|----------|-------|
| 11-0042  | 130°  |
| 11-0043  | 135°  |
| 11-0044  | 140°  |
| 11-0045  | 145°  |
| 11-0046  | 150°  |



## Lag Screw Tap

Cat. No. 7111-0014



## CLASSIC Insertion Wrench

Cat. No. 11-0054



## AMBI Insertion Wrench

Cat. No. 11-0022



## Replacement Retaining Rod for AMBI and CLASSIC Insertion Wrenches

Cat. No. 7111-0024



## AMBI CLASSIC Centering Sleeve

Cat. No. 11-0029







Plate Tamper

Cat. No. 11-0020



Insertion/Removal Wrench

Cat. No. 11-5061



Replacement Retaining Rod for  
Insertion/Removal Wrench

Cat. No. 7111-5062




Cannulated Barrel Guide

Cat. No. 7111-0060



Insertion/Removal Wrench  
Centering Sleeve

Cat. No. 7111-0030



Cannulated Plate Tamper

Cat. No. 11-0903



Plate Clamp

Cat. No. 21-0204



Combination Drill Guide, 3.5mm

Cat. No. 11-0075



Pin Guide, 2.4mm

Cat. No. 7111-0105



Guide Pin, 2.4mm

Cat. No. 41-0236



Twist Drill, 3.5mm

Cat. No. 7111-0045



Bone Screw Tap for 4.5mm  
Self-Tapping Screws

Cat. No. 11-0077



Bone Screw Tap for 4.5mm  
Nonself-Tapping Screws

Cat. No. 7111-0070



Bone Screw Length Gauge

Cat. No. 41-3500



Screw Pickup

Cat. No. 7111-5085



Self-Holding Hex Screwdriver

Cat. No. 7111-0026



Hex Screwdriver

Cat. No. 11-5035



Cannulated Plate Tamper

Cat. No. 7111-0038







## Lag Screw Trephine

Cat. No. 11-0926



## Long Slot Drill Guide

Cat. No. 11-5043



## Drill Sleeve

3.5 mm

Cat. No. 7111-0123



## Twist Drill

4.5 mm

Cat. No. 7111-0027

## Obturator

Cat. No. 11-6500



## Supracondylar Pin Guides


| Cat. No. | Angle |
|----------|-------|
| 11-0018  | 90°   |
| 11-0019  | 95°   |

## CLASSIC and AMBI Template

Cat. No. 7118-0454

(Not Shown)







Sterilization Tray  
Cat. No. 7111-5090



Standard Tray Insert  
Cat. No. 7111-5070



AMBI CLASSIC  
Instrument Tray  
Cat. No. 7111-5091



Screw Caddy  
Cat. No. 7111-5097



Bone Screw Caddy (optional)  
Cat. No. 7111-0137



Trial Tray (optional)  
Cat. No. 7111-5032













**Orthopaedics**

Smith & Nephew, Inc.  
1450 Brooks Road  
Memphis, TN 38116  
USA

[www.smith-nephew.com](http://www.smith-nephew.com)

Telephone 901-396-2121  
Information 1-800-821-5700  
Orders/Inquiries 1-800-238-7538