

Anatomical Shoulder™ Fracture

Surgical Technique



Anatomical Meets Fractures



Disclaimer

This document is intended exclusively for experts in the field, physicians in particular, and it is not intended for laypersons.

Information on the products and procedures contained in this document is of a general nature and does not represent medical advice or recommendations. Because this information does not constitute any diagnostic or therapeutic statement with regard to any individual medical case, examination and advising of the respective patient are absolutely necessary and are not replaced by this document in whole or in part.

Information contained in this document was gathered and compiled by medical experts and qualified Zimmer employees to the best of their knowledge. The greatest care was taken to ensure the accuracy and approachability of the information used and presented. Zimmer does not assume any liability, however, for the accuracy, completeness or quality of the information Zimmer is not liable for tangible or intangible losses that may be caused by the use of this information.

Surgical Technique Anatomical Shoulder Fracture

Developed in conjunction with

C. Gerber, MD, Prof. Zurich, Switzerland

J. JP. Warner, MD Boston, USA

L. Higgins, MD Boston, USA

A. Rüter, MD, Prof. Augsburg, Germany

Ph. Hardy, MD, Prof. Paris, France

Table of Contents

Foreword	4
Indications	4
Description of the Implants	5
Overview of the Instruments	8
Preoperative Planning	9
Surgical Technique	10
Patient Positioning and Surgical Approach	10
Delto-Pectoral Approach	10
Identification of the Lesser and Greater Tuberosities	10
Humeral Head Excision	11
Humeral Shaft Preparation	11
Sizing Convention and Consideration	12
Optional Retroversion Adjustment Technique	12
Optional Height Adjustment Technique	12
Mobilizing the Tuberosities	13
Assembling the Anatomical Shoulder Fracture Implant	16
Cementing the Prosthesis	17
Reattach Tuberosities	18
Closure	19
Ordering Information	21
Implants Instruments	21 22

Foreword

Shoulder replacement has evolved to achieve as a biomechanical and anatomical reconstruction of the shoulder.

Precise anatomical reconstruction using the *Anatomical Shoulder* System allows the surgeon to restore the geometry of the normal joint, thus ensuring good motion and pain relief, as well as durability of the reconstruction.

The Anatomical Shoulder Fracture System is specialized for anatomical reconstruction in case of proximal three or four part fractures.

The goal of a hemiarthroplasty for fracture is to replace the humeral head with a prosthetic component and to restore rotator cuff function by reconstructing the tuberosities to both, the shaft and the prosthesis.

Peter Wendt, Director Development Extremities

Markus Rauscher, Director Brand Management & Marketing Extremities

Indications

The *Anatomical Shoulder* Fracture System is indicated for:

- 4-part fracture of the proximal humerus
- 3-part or head split fractures

The goal of a hemiarthroplasty for fracture is to replace the humeral head with a prosthetic component and to restore rotator cuff function by reconstructing the tuberosity to both, the shaft and the prosthesis.

Description of the Implants

Specifically designed for anatomical reconstruction of articular proximal humerus fractures.

Anatomical Shoulder Fracture Head

Anatomical Shoulder Fracture Screw

Anatomical Shoulder Fracture Base-Plate

Anatomical Shoulder Fracture Stem

(FEI)

Anatomical Shoulder Fracture Base-Plate, Head and Screw

6 differ	ent versions	delivered	sterile in	one box

Right	40R, 44R, 48R
Left	40L, 44L, 48L



12 different versions delivered sterile in one box

Stem lengths	Short	Long	
	130 mm	170 mm	200 mm
Stem sizes	7, 8, 9, 10, 11, 12, 13, 14	7	9, 11, 13



Special Tuberosity Groove

Anatomical Design

- Anatomically designed proximal partThe proximal volume of the
- prosthesis is optimized to restore the normal given humeral anatomy
- Round shaped heads Special tuberosity groove
- Special tuberosity grooveBone fits anatomically below the head
- Prevention of impingement and tears



Special Fx Spikes

Special Fx Spikes

- For a fix and stable anchoring of the tuberosities to the stem
- Increase primary stability



Special Medial Clearance in Order to Preserve Remaining Humeral Bone

Effective Medial Base-Plate Design

• The base-plate has a special medial clearance in order to preserve remaining humeral bone



Extra Rough Blasted Proximal Surface

Surface

- Extrarough-blasted proximal surface with Fx spikes specially developed for proximal humeral bone
- Smooth-blasted distal surface
- Material: Optimized material composition, Ti Alloy, CoCr, with historical long-time experience
- Allows for cemented or noncemented use



Special Fx Suture Holes

Special Fx Suture Holes

- Due to the position of the special Fx suture holes, an anatomical repositioning of the greater and lesser tuberosity below the head, back to the original anatomy is possible.
- Therefore reconstruction failures are minimized.
- Fx suture holes optimized for secure and safety fixation.
- High-end fluid ground and rounding Fx suture holes



- To provide final stable fixation of the tuberosities to the Anatomical Shoulder Fracture Stem.
- One medial, two additional lateral



Special Fx Suture Holes



Safety Distal Anchoring

Safety Distal Anchoring

- 3 grooves for rotational stability, press-fit or cemented
- The two lateral grooves are for defined orientation
- Tapered stem design, for stem removal and forces are distributed homogeneously



Convertible to an Anatomical Shoulder Inverse/Reverse System, without Stem Removal

Such revision might be necessary in case of irreparable rotator cuff tear. This will greatly simplify and shorten revision surgery since the need to remove a well-fixed stem is eliminated.



Right and Left Option

Right and Left Option

- In respect to the shoulder anatomy, physiology, and the greater and lesser tuberosities
- Continuous adjustment of retroversion



Anatomical Shoulder Fracture Tray

The preparation and implantation of the *Anatomical Shoulder* Fracture System should be carried out in a standardized manner. The special set of instruments has been logically developed. The required instruments have been limited to a minimum. The correct use and handling of these special devices are a requirement of the success of the surgery. To convert an *Anatomical Shoulder* Fracture System into an *Anatomical Shoulder* Inverse/Reverse System, without the need of stem removal, further instruments are needed:

- Anatomical Shoulder Inverse/ Reverse Tray
- Anatomical Shoulder Instrument Tray I and II
- Anatomical Shoulder Glenoid Tray

See Surgical Technique Anatomical Shoulder Inverse/Reverse (Lit. No. 06.01276.012).



Lit. No. 06.01276.012



 Anatomical Shoulder Fracture Long Stem Lit. No. 06.01372.000

The following radiographic images of the shoulder joint are desired for preoperative planning:

- Full-size true anterior-posterior view with neutral rotation (0°), centered on the articular cavity
- Axial view
- Y view
- CT scan

An initial assessment is made of the bone in the superior and inferior aspects of the shoulder, using radiographic and CT imaging in order to determine the suitability of the patient's available bone stock for implant insertion.

Preoperative planning is also carried out, using AP and lateral shoulder radiographs of known magnification, and the available templates to confirm the size and alignment of the implant.

Surgical Technique

Patient Positioning and Surgical Approach

The patient should be placed in a "beach chair" position on the edge of the operating table (Fig. 1).

The arm must be freely movable and it must be possible to extend it fully. An armrest is optional.



Delto-Pectoral Approach

Make a skin incision in a straight line starting from the lateral edge of the coracoid as far as the insertion of the deltoid muscle. Seek the cephalic vein between the deltoid muscle and the pectoralis major muscle. Make the approach medial to the vein, to open the deltopectoral groove.

The coracoid process is identified. The clavi-pectoral fascia is incised at the external border of the coracobrachialis. The axillary nerve is then identified before identification of the subscapularis.

In fracture cases, it is especially important to identify and protect the musculocutaneous and the axillary nerves.

Identification of the Lesser and Greater Tuberosities

The glenohumeral joint is exposed by extending the fracture line between the tuberosities, incising the rotator interval over the long head of the biceps tendon. The biceps tendon is an excellent landmark to identify the interval between the lesser, and the greater tuberosity. Even if the biceps tendon has been ruptured, place the scissors in the bicipital groove and use them to open the interval between the subscapularis and the supraspinatus tendon. Next, free up the lesser tuberosity from the underlying humeral head and soft tissues. Now, in a similar manner, carefully identify and free up the greater tuberosity.

The greater and lesser tuberosity fragments must be sufficiently freed up so that they can be easily repaired around the *Anatomical Shoulder* Fracture and to each other at the time of closure.

Humeral Head Excision

With the tuberosities retracted out of the way, use a clamp to retrieve the humeral head.

Now compare its dimension with one of the 3 head sizes of the *Anatomical Shoulder* Fracture System.

If the humeral head is between available prosthetic head sized, select the smaller of the two. The most common mistake is to use a too large humeral head (Fig. 8).

Humeral Shaft Preparation

Insert the locking spring for rasp into the lateral slot of the rasp, for stabilizing height positioning during preparation (Fig. **3**).

Now attach the Handle for rasp to the rasp (Fig. 4).

Manually rasp the humeral canal using progressively larger rasps in 1 mm increments until slight resistance is felt from cortical contact in the canal.

Ream to the appropriate depth for the selected stem lengths. The depth corresponds to the implant length to be used. If a long *Anatomical Shoulder* Fracture Stem is required, connect a rasp extension on the distal end of the rasp.

Rasp extensions are available for *Anatomical Shoulder* Fracture Stems sizes 7, 9, 11 and 13.



Optional Retroversion Adjustment Technique

Insert the alignment rod into the appropriate retroversion hole on the handle for rasp. Use the right or left hole for the corresponding shoulder side and the preferred hole for orientation to the forearm or to the condyles (Fig. 5).

Optional Height Adjustment Technique

The outer-shaped *Anatomical Shoulder* Fracture Heads are laser marked on the handle for rasp for height orientation during rasp procedure. The correct rasp deepness is reached, if you feel that the laser-marked head is in right high position (Fig. 6).

Attach the *Anatomical Shoulder* Fx Ruler to the Handel for Rasp for your height adjustment control. Now, use the pectoralis for height orientation. On the *Anatomical Shoulder* Fx Ruler you will find a laser-marked area of the upper border of the pectoralis major tendon.

Now control if the laser-marked head on the handle is placed in the right height and the laser-marked area of the pectoralis corresponds to the **upper border of the pectoralis.**

Sizing Convention and Consideration

AS Fx rasp size	Rasp extension	Cemented implant size and length	Press-fit, implant size and length
7р	yes	-	7-130, 7-170
8p	no	-	8-130
9p/7c	yes	7-130, 7-170	9-130, 9-200
10p/8c	no	8-130	10-130
11p/9c	yes	9-130, 9-200	11-130, 11-200
12p/10c	no	10-130	12-130
13p/11c	yes	11-130, 11-200	13-130, 13-200
14p/12c	no	12-130	14-130
13c	yes	13-130, 13-200	-
14c	no	14-130	-

Fig. 5

ł

p = press-fit (uncemented) c = cemented



Disconnect the handle for rasp and insert the screw for rasp with the help of the screwdriver.

The screw will press against the spring and the spring against the bone. Now, the rasp is fixed in the appropriated position, in height and retroversion, in the humeral shaft (Fig. 7).

If the humeral head is between available prosthetic head sized, select the smaller of the two. The most common mistake is to use a too large humeral head (Fig. 8).

Use the right or left humeral trial head component for the corresponding shoulder side.

Then attach the selected trail head to the rasp, seated in the humeral shaft. With the screw for humeral trial head and the screwdriver prepare the stable fixation of the trial components (Fig. 9).



Fig. 8

Mobilizing the Tuberosities

When the proper height and torsion of the trial prosthesis has been determined, mobilize the tuberosities so that they can be approximated around the prosthesis, to one another and to the humeral shaft. Due to the position of the special Fx suture holes, an anatomical repositioning of the greater and lesser tuberosity below the head, back to the original anatomy is possible (Fig. 9).



The primary goal of tuberosity reattachment is to get as much contact with the stem and the proximal humeral shaft while rebuilding into the anatomical position.

The initial reduction of the greater tuberosity enables both the height and the retroversion to be tested. The greater tuberosity is placed on the diaphysis and the prosthesis, use the special Fx suture hole for the greater tuberosity and place it in the tuberosity groove under the round-shaped head (Fig. 10).

Now, the height is tested to ensure it is correct:

- The tension of the supraspinatus and the long head of biceps which must arch over the *Anatomical Shoulder* Fracture Head, and the height of the acromio-humeral space.
- The top of the greater tuberosity must be located below the upper pole of the *Anatomical Shoulder* Fracture Head.
- There must be no diastasis or overlap between the greater tuberosity and the humeral diaphysis.

Now test the retroversion:

• Arm in neutral position – the *Anatomical Shoulder* Fracture Head must face the glenoid.

When the version and height of the *Anatomical Shoulder* Fracture Trails (rasp and trail head) are set, insert the *Anatomical Shoulder* Fx Ruler into the proximal-lateral suture hole of the head component and mark the position next to the associate laser mark on the *Anatomical Shoulder* Fx Ruler.

After reducing the joint, perform a final range of motion assessment (Fig. 10).



Now, remove the ruler, unscrew the *Anatomical Shoulder* Fracture Head Trial, remove the screw for rasp, attach the handle for rasp to the rasp and remove the rasp from the humeral shaft.

Clean the fracture site at the shaft edges and place drill holes through the shaft; two lateral and two medial to the biceps groove. Place sutures through the shaft drill holes. These will be greater tuberosity and lesser tuberosity vertical sutures that will go up around the top of the bone segments through the rotator cuff-bone junction (Fig. **11**).



Assembling the Anatomical Shoulder Fracture Implant

Humeral stem implant size is selected based upon technique and fixation desired. For example, choose for rasp size 11p/9c the implant stem 9 cemented. If a press-fit is desired, choose for rasp size 11p/9c the implant stem 11 (refer to Sizing Convention and Consideration section for additional information, page 12). Humeral head implant size is the same size and version (left or right) as the trial head chosen previously.

Connect and assemble the elected Anatomical Shoulder Fracture Base-Plate to the Anatomical Shoulder Fracture Stem with the locking screw (Fig. 12).

Note: Due to the ability to convert from anatomical to inverse/reverse, there has to be a gap between stem and base-plate.

First, connect the distal support to the Fx adapter (A).

Place the assembled stem into the special stem holder (Fig. 13) and close the adpter arm (B).

Connect finally with the help of the torque wrench, stem with baseplate (Fig. 14).



Complete the *Anatomical Shoulder* Fracture Implant by impacting the *Anatomical Shoulder* Fracture Head to the stem (Fig. **15**).

Attach the *Anatomical Shoulder* Fx Ruler onto the *Anatomical Shoulder* Fracture Implant, by inserting the peg of the ruler into the superior hole of the base-plate.



Cementing the Prosthesis

Thoroughly irrigate the medullary canal to remove blood and other debris. Insert a cement blocker at the appropriate depth in the medullary canal. If possible use high-viscosity cement mixed under vacuum and insert it with a cement gun.

Insert the *Anatomical Shoulder* Fracture Implant into the humeral canal to the same level of the *Anatomical Shoulder* Fx Ruler relative to the mark, marked earlier when the rasp and the trial head were used (Fig. **16). Use** the head impactor for finally impaction.

If cemented, make certain that there is no excess cement extruding from the canal proximally above the humeral stem into the fracture site. This will interfere with the potential for bony union between the tuberosities, stem, and the diaphyseal fragment. Use a curette to remove any excess cement. It is important to keep the sutures separated to avoid confusion in tying the proper sutures (Fig. 16).

Reattach Tuberosities

Fixation of the tuberosities is critical to the success of the procedure. Basic principles in fracture repair should be followed to provide stable fixation of the tuberosities into the stem. The following description provides guidelines for using of the suture holes in the stem to provide proper fixation. Suture pattern and method can be modified based on the condition of the fracture. The primary goal of tuberosity reattachment is to get as much contact with the stem and the proximal humeral shaft while rebuilding into the anatomical position.

A suture should be placed in the special greater tuberosity suture hole (green), a second suture in the special lesser tuberosity suture hole (orange) (Fig. 16). These sutures will initially be used to position the tuberosities to the shaft in a cerclage fashion (Fig. 17).

The posterior end of the suture, passed to the greater tuberosity suture channel (green), is passed at the junction between posterior end of the supraspinatus tendon and greater tuberosity, the anterior part is passed inside out at the junction greater tuberosity and anterior border of the supraspinatus tendon. This suture is then tied and reduces the greater tuberosity anatomically in an anatomic fashion.

The cerclage sutures will be passed through the subscapularis tendon, at its insertion, wrap around the lesser and greater tuberosities and pass through the infraspinatus and teres minor at the tendon insertions. These sutures will be tightened and tied off first.

The sutures placed in the humeral shaft lateral to the biceps groove, will be passed through the supraspinatus tendon at its insertion and used to bring the distal edge of the greater tuberosity back down to the shaft.



The sutures placed in the humeral shaft medial to the biceps groove, will be passed through the subscapularis at the tendon insertion and used to bring the distal edge of the lesser tuberosity back down to the shaft. These vertical sutures will be tightened and secured after the cerclage sutures are tied off.

The sutures placed in the middle of the *Anatomical Shoulder* Fracture Implant are used to further reduce or compress the fragments against the prosthesis, if necessary.

A suture from each hole will be passed posteriorly through Infraspinatus and Teres minor insertions, respectively. The suture exiting anteriorly will pass around the greater tuberosity fragment and be tied down onto the greater tuberosity.

If necessary, a second suture from each hole will be passed posteriorly around the stem and the medial hole through the subscapularis at its insertion. The suture end exiting anteriorly will be wrapped around the lesser tuberosity and tied down against the lesser tuberosity (Fig. 18 and 19).

Remove and discard any unused sutures. Close the rotator interval from the edge of the supraspinatus to the upper edge of the subscapularis tendon.

Check stability and range of motion. If necessary, place bone graft from the humeral head in and around the tuberosity shaft interface.

Closure

Close the subcutaneous layers, and then the skin.





Anatomical Shoulder[™] Fracture Implants



Anatomical Shoulder™ Fracture Stems Cemented/Uncemented



Anatomical Shoulder™ Fracture Long Stems Cemented/Uncemented

Size	L [mm]	REF
7	170	01.04217.072
9	200	01.04217.092
11	200	01.04217.112
13	200	01.04217.132



Anatomical Shoulder™ Fracture Heads

Anatomical Shoulder™ Fracture Head Protasul®-21WF ISO 5832-12

S

Anatomical Shoulder™ Fracture Screw Protasul®-21WF ISO 5832-12



Anatomical Shoulder™ Fracture Base-Plate Protasul®-100 ISO 5832-11

STERILE R

Size (S)	REF	
40	Left	01.04227.400
40	Right	01.04227.405
44	Left	01.04227.440
44	Right	01.04227.445
48	Left	01.04227.480
48	Right	01.04227.485

Instruments for Anatomical Shoulder[™] Fracture

 Article
 REF

 Anatomical Shoulder™

 Fracture Instrument Tray

 (complete)
 ZS01.04237.000

Tray Cover 01.00029.031

Anatomical Shoulder™ Fx Tray (empty) 01.04237.010

Anatomical Shoulder™ Fx Insert I for Tray (empty) 01.04237.020

Anatomical Shoulder™ Fx Humeral Trial Heads

Left 48	
Right 48	
Left 44	
Right 44	
Left 40	
Right 40	

01.04237.480 01.04237.485 01.04237.440 01.04237.440 01.04237.400 01.04237.405

Anatomical Shoulder[™] Fx Adapter 01.04237.600



Article REF **Torque Wrench for Humeral Head** 72.11.20-06

Anatomical Shoulder[™] Fx Ruler 01.04237.700

Anatomical Shoulder™ Fx Screw for Humeral Trial Head 01.04237.500

Nut for Torque Wrench, 4,5 mm 72.11.20-09

Impactor for Humeral Stem 72.01.00-01

Hexagonal Wrench, 5 mm 5331

Anatomical Shoulder™ Fx Humeral Stem Setting Instrument 01.04237.610

Hexagonal Screwdriver, 2.5 mm 109.02.020

Article	REF
Anatomical Shoulder ^T	™ Fx Insert II
for Tray (empty)	01.04237.030

Anatomical Shoulder[™] Fx Rasps (Set 1)

01.04237.075
01.04237.095
01.04237.115
01.04237.135
01.04237.155

Anatomical Shoulder[™] Fx Rasps (Set 2)

Size 6.5 [7p] Size 8.5 [9p/7c] Size 10.5 [11p/9c] Size 12.5 [13p/11c] 01.04237.125 Size 14.5 [13c]

01.04237.065 01.04237.085 01.04237.105 01.04237.145

Anatomical Shoulder[™] Fx 01.04237.180

Springs for Rasp



Article REF Anatomical Shoulder™ Fx Handle for Rasp 01.04237.200

Anatomical Shoulder[™] Fx **Rasp Extensions**

Size 6.5 [7p]	01.04237.066
Size 8.5 [9p/7c]	01.04237.086
Size 10.5 [11p/9c]	01.04237.106
Size 12.5 [13p/11c]	01.04237.126
Size 14.5 [13c]	01.04237.146

Anatomical Shoulder™ Fx Screws for Rasp 01.04237.190

Anatomical Shoulder™ Fx Alignment Rod 01.04237.310



Contact your Zimmer representative or visit us at www.zimmer.com



