



# ANATOMICALLY CONTOURED

With Low-Profile Fixation

# DISTAL VOLAR RADIAL RADIUS PLATES

Available with 2, 3, 4, 6, 8, 10, and 12 holes in both right and left configurations (Use 2.7mm, 3.5mm and 4.0mm screws).

# DISTAL VOLAR ULNAR RADIUS PLATES

Available with 2, 3, 5, and 7 holes in both right and left configurations (Use 2.7mm, 3.5mm and 4.0mm screws).





### **DISTAL VOLAR RADIAL PLATES**

## **Product Information**

Distal Volar Radial Radius Plates: Stainless Steel (Volar-plate orientation towards the Radius)

Distal Volar Ulnar Radius Plates: Stainless Steel (Volar-plate orientation towards the Ulna)


### For use with:

- Periarticular 3.5mm Self-Tapping Cortical Screws with 2.7mm head
- 2.7mm Mini fragment Cortical Screws
- Forte® Distal Radial Plate System 4.0mm partially or fully threaded Cancellous Screws



### REDUCED PROFILE SCREWS

The head size of the Periarticular screws has been reduced significantly to create a nearly flush profile when used with a plate or independently.



**DISTAL VOLAR** 

DISTAL VOLAR RADIAL PLATES	
DISTAL VOLAR RADIAL RADIUS PLA	TE
Plate	
Distal Volar Radial Radius Plate	

DIGHTE VOLIK KILDINE KILDIOG I LITEG				
Plate	Left	Right	Plate Holes	Length
Distal Volar Radial Radius Plate	2348-014-02	2348-013-02	2	32mm
Distal Volar Radial Radius Plate	2348-014-03	2348-013-03	3	43mm
Distal Volar Radial Radius Plate	2348-014-04	2348-013-04	4	54mm
Distal Volar Radial Radius Plate	2348-014-06	2348-013-06	6	74mm
Distal Volar Radial Radius Plate	2348-014-08	2348-013-08	8	95mm
Distal Volar Radial Radius Plate	2348-014-10	2348-013-10	10	116mm
Distal Volar Radial Radius Plate	2348-014-12	2348-013-12	12	137mm

# DISTAL VOLAR LILNAR RADIUS PLATES

Plate	Left	Right	Plate Holes	Length
Distal Volar Ulnar Radius Plate	2348-018-02	2348-017-02	2	33mm
Distal Volar Ulnar Radius Plate	2348-018-03	2348-017-03	3	44mm
Distal Volar Ulnar Radius Plate	2348-018-05	2348-017-05	5	64mm
Distal Volar Ulnar Radius Plate	2348-018-07	2348-017-07	7	85mm

	DISTAI	VOLAR
ULNAR	RADIU	s plate

BONE SCREW OPTIONS			
Screw Type	Cat No.	Length Options (mm)	
2.7mm Self-Tapping Cortical Screws	4827-mm-01	6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65, 70	
3.5mm Self-Tapping Cortical Screws with 2.7mm head	2348-mm-35	8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110	
Forte 4.0mm Partially Threaded Cancellous Screws	2480-mm-40	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34 36, 38, 40	
Forte 4.0mm Fully Threaded Cancellous Screws	2480-mm-41	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40	

	/ 11016	7 Halo
Ö	ט דוסומ	л LOO
0	٦	5 [
	2 HOIE	3 Holo

SCREW ANGULATION		
Screw Type	Shaft Angulation (Dua Longitudinal	al Compression Slot) Transverse
2.7mm Cortical Screw	41°	19°
3.5mm Cortical Screw (w/ 2.7mm)	37°	15°

#### For use with the following sets:

- Periarticular 2.7mm/3.5mm/4.0mm Screw and Instrument Set #2348-000-00
- · Mini and Small Fragment Instrument and Implant sets
- Distal Radial Plates Case #2348-060-00
- Distal Radial Plates Standard Set #2348-000-03
- Distal Radial Plates Expanded Set #2348-000-04



For more information about the Periarticular Plating System, contact your Zimmer representative, or visit us at www.zimmer.com.