



ANATOMICALLY CONTOURED

With Low-Profile Fixation

DISTAL LATERAL FEMORAL PLATES

4.5mm Plate: available with 6, 8, 10, 12, 14, 16, and 18 holes in both right and left configurations (Use 4.5mm and 6.5mm screws).





DISTAL LATERAL FEMORAL PLATE

Product Information

DISTAL LATERAL FEMORAL PLATES

Distal Lateral Femoral Plate

Distal Lateral Femoral Plate

Distal Lateral Femoral Plate

Distal Lateral Femoral Plates: Stainless Steel

For use with:

Plate

 Periarticular 4.5mm Self-Tapping Cortical Screws, and 6.5mm Cancellous Screws



REDUCED PROFILE SCREWS

Plate Holes

6

8

10

The head size of the Periarticular screws has been reduced significantly to create a nearly flush profile when used with a plate or independently.

Length

146mm

178mm

210mm

18 Hole—
—16 Hole—
—14 Hole —
— 12 Hole —
—10 Hole —
— 8 Hole —
6 Hole

Distal Lateral Femoral Plate	2347-02-12	2347-01-12	12	242mm
Distal Lateral Femoral Plate	2347-02-14	2347-01-14	14	273mm
Distal Lateral Femoral Plate	2347-02-16	2347-01-16	16	305mm
Distal Lateral Femoral Plate	2347-02-18	2347-01-18	18	337mm
	1			
SCREW ANGULATION				

Left

2347-02-06

2347-02-08

2347-02-10

Right

2347-01-06

2347-01-08

2347-01-10

SCREW ANGULATION		
Screw Type	Shaft Angulation (Dual Compression Slot) Longitudinal Transverse	
4.5mm Periarticular Cortical Screw	62°	12°

Screw Type	Cat No.	Length Options (mm)
3.5mm Self-Tapping Cortical Screws	2347-20-mm 2347-120-mm	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150
4.0mm Cancellous Screws	2347-21-mm (partially threaded) 2347-22-mm (fully threaded)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 35, 40, 45, 50, 55, 60 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 45, 50, 55, 60
4.5mm Self-Tapping Cortical Screws	2347-23-mm 2347-123-mm	26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150
6.5mm Cancellous Screws	2347-24-mm 2347-124-mm (with 16mm thread) 2347-25-mm 2347-125-mm (with 32mm thread)	30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150
	2347-26-mm 2347-126-mm (fully threaded)	25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 95, 100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150

For use with the following sets:

- Periarticular Instrument Set #2347-000-04
- Periarticular Screw Tray #2348-45 (For 3.5mm, 4.0mm, 4.5mm & 6.5mm screws)
- Distal Femoral Plates Case #2348-080-00
- Distal Femoral Standard Plate Set #2347-000-12 (Includes 4 and 6 Holes, left and right Medial and 6, 8, and 10 Holes left and right Lateral)
- Distal Femoral Expanded Plate Set #2347-000-13 (All Distal Femoral Plates)





For more information about the Periarticular Plating System, contact your Zimmer representative, or visit us at www.zimmer.com.