

PERIARTICULAR Plating system

ANATOMICALLY CONTOURED

With Low-Profile Fixation

PROXIMAL LATERAL TIBIAL PLATES

3.5mm Plate: available with 4, 6, 8, 10, 12, 14, and 16 holes in both right and left configurations (Use 3.5mm and 4.0mm screws).

4.5mm Plate: available with 2, 4, 6, 8, 10, 12, and 14 holes in both right and left configurations (Use 4.5mm and 6.5mm screws).

3.5mm Plate





PROXIMAL LATERAL TIBIAL PLATE

Product Information

Proximal Lateral Tibial Plates: Stainless Steel

For use with:

- Periarticular 3.5mm and 4.5mm Self-Tapping Cortical Screws, 4.0mm and 6.5mm Cancellous Screws
- Small Fragment Set for 3.5mm plates and Basic Fragment Set for 4.5mm plates



REDUCED PROFILE SCREWS

The head size of the Periarticular screws has been reduced significantly to create a nearly flush profile when used with a plate or independently.

3.5mm PROX	IMAL LATERAL	TIBIAL PLATE		4.5mm PROXIMAL LATERAL TIBIAL PLATE			
Left	Right	Plate Holes	Length	Left	Right	Plate Holes	Length
_	_	_	_	2347-06-02	2347-05-02	2	42mm
2347-04-04	2347-03-04	4	59mm	2347-06-04	2347-05-04	4	73mm
2347-04-06	2347-03-06	6	83mm	2347-06-06	2347-05-06	6	105mm
2347-04-08	2347-03-08	8	107mm	2347-06-08	2347-05-08	8	137mm
2347-04-10	2347-03-10	10	131mm	2347-06-10	2347-05-10	10	169mm
2347-04-12	2347-03-12	12	154mm	2347-06-12	2347-05-12	12	201mm
2347-04-14	2347-03-14	14	182mm	2347-06-14	2347-05-14	14	232mm
2347-04-16	2347-03-16	16	202mm	_	_	_	_

Screw Type	Cat No.	Length Options (mm)
3.5mm Self-Tapping Cortical Screws	2347-20-mm	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95
	2347-120-mm	100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150
4.0mm Cancellous Screws	2347-21-mm (partially threaded)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 35, 40, 45, 50, 55, 60
	2347-22-mm (fully threaded)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 45, 50, 55, 60
4.5mm Self-Tapping Cortical Screws	2347-23-mm	26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95
	2347-123-mm	100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150
6.5mm Cancellous Screws	2347-24-mm (with 16mm thread)	30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95
	2347-124-mm (with 16mm thread)	100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150
	2347-25-mm (with 32mm thread)	45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95
	2347-125-mm (with 32mm thread)	100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150
	2347-26-mm (fully threaded)	25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 95
	2347-126-mm (fully threaded)	100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150


For use with the following instrument sets:

- Periarticular Instrument Case #2348-50
- Periarticular Screw Case #2348-45 (3.5mm, 4.0mm, 4.5mm and 6.5mm screws)



Contact your Zimmer representative, or visit us at www.zimmer.com.

3.5mm Plate



4.5mm Plate

