



ANATOMICALLY CONTOURED

With Low-Profile Fixation

LATERAL CALCANEAL PLATE

Allows use of 2.7mm screws, 3.5mm screws with 2.7mm head, and 4.0mm *Forte* screws both fully threaded and partially threaded (Available in left and right plates).





LATERAL CALCANEAL PLATE

Product Information

Lateral Calcaneal Plates: Stainless Steel


For use with:

- Periarticular 3.5mm Self-Tapping Cortical Screws with 2.7mm head
- 2.7mm Mini Fragment Cortical Screws
- Forte® Distal Radial Plate System 4.0mm partially or fully threaded Cancellous Screws



REDUCED PROFILE SCREWS

The head size of the Periarticular screws has been reduced significantly to create a nearly flush profile when used with a plate or independently.



BONE SCREW OPTIONS		
Screw Type	Cat No.	Length Options (mm)
2.7mm Self-Tapping Cortical Screws	4827-mm-01	6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65, 70
3.5mm Self-Tapping Cortical Screws with 2.7mm head	2348-mm-35	8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110
Forte 4.0mm Partially Threaded Cancellous Screws	2480-mm-40	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40
Forte 4.0mm Fully Threaded Cancellous Screws	2480-mm-41	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40

For use with the following instrument sets:

- Periarticular 2.7mm/3.5mm/4.0mm Screw and Instrument Set #2348-000-00
- Mini and Small Fragment Instrument Sets



For more information about the Periarticular Plating System, contact your Zimmer representative, or visit us at www.zimmer.com.