

ANATOMICALLY CONTOURED

**With Low-Profile
Fixation**

LATERAL DISTAL TIBIAL PLATES

Available with 4, 6, 8, 10, 12, 14, 16, 18, and 20 holes in both right and left configurations (Use 3.5mm and 4.0mm screws).

MEDIAL DISTAL TIBIAL PLATES

Available with 4, 6, 8, 10, 12, 14, 16, 18, and 20 holes in both right and left configurations (Use 3.5mm and 4.0mm screws).

The periarticular screws feature a reduced profile head which is designed to create a nearly flush profile when used with a plate or independently

Compression screw holes will accommodate periarticular screws or conventional stainless steel cortical screws

Reduced plate thickness in periarticular regions helps to minimize soft tissue irritation

Parallel placement of screws to joint line is possible while allowing angulation if necessary

Lateral Distal Tibial Plates

Medial Distal Tibial Plates



Preoperative



Postoperative



Preoperative



Postoperative



LATERAL & MEDIAL DISTAL TIBIAL PLATES

Product Information

Lateral & Medial Distal Tibial Plates: Stainless Steel

For use with:

- Periarticular 3.5mm Self-Tapping Cortical screws, 4.0mm Cancellous Screws
- Small Fragment Set for 3.5mm plates



REDUCED PROFILE SCREWS

The head size of the Periarticular screws has been reduced significantly to create a nearly flush profile when used with a plate or independently.

LATERAL DISTAL TIBIAL PLATE				MEDIAL DISTAL TIBIAL PLATE			
Left	Right	Plate Holes	Length	Left	Right	Plate Holes	Length
2347-10-04	2347-09-04	4	61mm	2347-08-04	2347-07-04	4	96mm
2347-10-06	2347-09-06	6	85mm	2347-08-06	2347-07-06	6	120mm
2347-10-08	2347-09-08	8	109mm	2347-08-08	2347-07-08	8	144mm
2347-10-10	2347-09-10	10	133mm	2347-08-10	2347-07-10	10	168mm
2347-10-12	2347-09-12	12	157mm	2347-08-12	2347-07-12	12	192mm
2347-10-14	2347-09-14	14	181mm	2347-08-14	2347-07-14	14	216mm
2347-10-16	2347-09-16	16	205mm	2347-08-16	2347-07-16	16	240mm
2347-10-18	2347-09-18	18	229mm	2347-08-18	2347-07-18	18	264mm
2347-10-20	2347-09-20	20	253mm	2347-08-20	2347-07-20	20	288mm

BONE SCREW OPTIONS

Screw Type	Cat. No.	Length Options (mm)
3.5mm Self-Tapping Cortical Screws	2347-20-mm	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95
	2347-120-mm	100, 105, 110, 115, 120, 125, 130, 135, 140, 145, 150
4.0mm Cancellous Screws	2347-21-mm (partially threaded)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 35, 40, 45, 50, 55, 60
	2347-22-mm (with full thread)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 45, 50, 55, 60

For use with the following sets:

- Periarticular Distal Tibial Plate Standard Set #2347-000-10 (Includes 4, 6, 8, and 10 hole plates)
- Periarticular Distal Tibial Expanded Set #2347-000-11 (Includes all Distal Tibial plates)
- Periarticular Distal Tibial Plate Case #2348-75
- Periarticular Instrument Tray #2347-000-04
- Periarticular Screw Tray #2348-45




zimmer
Confidence in your hands™



Contact your Zimmer representative, or visit us at www.zimmer.com.

Lateral Distal Tibial Plate



Medial Distal Tibial Plate

